EXA'lence™ 370

VINYL POLYETHER SILICONE (VPES™) IMPRESSION MATERIAL
For use only by a dental professional in the recommended indications.

Viscosity	Setting	Type of Mixing Tip
MEDIUM BODY (MONOPHASE)	Regular & Fast Set	Dynamic mixing tip (yellow)
HEAVY BODY	Regular & Fast Set	Dynamic mixing tip (yellow)

RECOMMENDED INDICATIONS and impression techniques

Indication	Impression Technique	Tray Type	Recommended Viscosities
Impression for facing, inlay, onlay, crown or bridge	Single mix single impression	Stock or Custom Tray	Medium Body
	Double mix single impression		Heavy Body or Medium Body + (Extra) Light Body
	Double mix double impression (Putty-wash impression)		Putty or Heavy Body + Medium Body or (Extra) Light Body
	Triple tray technique	Triple Function Tray	Heavy Body (Rigid) + (Extra) Light Body
Impression for partial denture	Single mix single impression	Custom Tray	Medium Body
Impression for full denture	Single mix single impression	Custom Tray	Medium Body or Light Body

CONTRAINDICATIONS

In rare cases, the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

PHYSICAL PROPERTIES (Average)

Test conditions: Temperature (74°F/23°C ± 4°F/2°C) Relative humidity (50 ± 5%)
(ISO 4823: 2000 E) (Dental elastomeric impression materials)

	Medium Body (Monophase) (Type 2: Medium Consistency)	Heavy Body (Type 1: High Consistency)
Total working time (minutes)	Regular: 2 min. Fast: 1 min	Regular: 2 min. Fast: 1 min
Setting time (minutes)	Regular: 5 min. Fast: 3.5 min.	Regular: 5 min. Fast: 3.5 min.
Minimum time in mouth (minutes)	Regular: 3 min. Fast: 1.5 min.	Regular: 3 min. Fast: 1.5 min.
Recovery from deformation (%)	Regular: 98.4% Fast: 97.0%	Regular: 98.1% Fast: 98.1%
Maximum strain in compression (%)	Regular: 2.9% Fast: 3.4%	Regular: 3.0% Fast: 3.2%
Linear dimensional change (%)	≤ - 0.10 %	≤ - 0.10 %
Detail reproduction	20µm	50µm
Compatibility with gypsum	20µm	50µm

PRECAUTIONS

When extruding EXA'lence, care should be taken to avoid mixing or contact with the following materials. These may delay or prevent setting: catalyst for condensation silicone impression materials, polysulfide impression materials, eugenol materials, sulfur, latex, oil, acrylates, glycerol, spray-type or ointment-type topical surface anesthesia (e.g. lidocaine).

DIRECTIONS FOR USE

CARTRIDGE LOADING AND DISPENSING

Note: The cartridges are heavy and dropping may result in damage to the cartridges rendering it unusable. The cartridge and mixing tip provided are compatible with a dynamic mixer.

Also refer to the manufacturer's instructions.

- Hold and support the cartridge with the extrusion nozzles facing upwards.
- Grasp the tab on the cap and bend in an upward 90° direction. Do not place excessive stress on tab.
- While holding the cartridge securely, grasp the tab with thumb and forefinger and remove the cap from the cartridge.
- Insert the cartridge into the mixer according to the manufacturer's instructions.
- Start the mixer according to the manufacturer's instructions. Bleed a small amount of material from the two nozzles of cartridge prior to each use. Wipe off the extruded material in a vertical direction to avoid cross contamination of the nozzles.
- Install a supplied mixing tip on the extrusion nozzles following the manufacturer's instructions. NOTE: Do not use other mixing tips. They will not fit properly and can cause material to leak from the nozzles.
- If the mixing tip does not seat easily, be sure the internal central hex of the mixing tip is properly aligned with the drive shaft of the mixer.
- Once the mixing tip is properly placed, place nozzle locking ring over the mixing tip and turn it 1/4 turn clockwise until ring is locked in place.

IMPRESSION STEP

Note:

Prior to use, leave the material to stand at room temperature. Cooled or refrigerated material may delay setting.

- Make sure the preparation is dry and free of debris during impression taking.
- Block out undercuts to aid impression removal.
- Apply a tray adhesive (GC Universal VPS Adhesive, etc.) onto the tray and leave to dry. Read manufacturer's instructions for use.
- Press the button and dispense the mixed base and catalyst into an intraoral syringe (Monophase type only) or suitable tray.
- Immediately after dispensing, move the plungers of the machine off the cartridge to prevent continued pressure on the material and extrusion of excess material from the tip.
- Inject the wash material onto the prepared teeth and completely cover the preparation. Blasting with air is not necessary. EXA'lence's hydrophilic property enables better self-wetting, allowing an excellent adaptation to the tooth.

Note: In case of using the double mix double impression or Putty-wash technique, make sure the preliminary impression is cleaned and thoroughly dried prior to taking the final impression.

- Seat the loaded tray in the mouth. See above table for recommended working times.
- Wait for the material to set. See above table for recommended intra oral setting times.
 - EXA'lence has unique chemistry that allows it to polymerize faster in the oral cavity than at room temperature. Therefore to determine setting, examine material in the mouth since material on bench will yield a false result. Material at room temperature will not set within acceptable time.
 - Higher temperatures shorten working time and lower temperatures increase working time.
- Remove impression and rinse under running tap water.

DISINFECTION

The impression can be disinfected with a glutaraldehyde solution or other appropriate disinfectants. Refer to the respective manufacturer's instructions for use.

MODEL FABRICATION

- After the impression is rinsed with water, wait for 60 minutes before pouring a model. If the impression is disinfected, wait for 60 minutes after disinfecting the impression before pouring the model.
- Pour a model within 14 days after impression taking. EXA'lence is compatible with type III and type IV stone (GC Fuji ROCK™ EP) and can be silver or copper plated.

CAUTION

- Avoid contact of base, catalyst, mixed material or adhesive with eyes. In case of contact, immediately flush with generous amount of water and seek medical attention.
- Avoid getting material on clothing as it is difficult to remove.

STORAGE

Store in a dry location at normal room temperature (59 - 77°F; 15 - 25°C) (Relative humidity: 50 +/-10%). Avoid direct sunlight.

Shelf life: 2 years from date of manufacture.

PACKAGES

- EXA'lence 370 MEDIUM BODY (MONOPHASE)
Regular and Fast Set
511g (370mL) x 1 cartridge, Mixing tip for 370mL cartridge (yellow) x 10 pieces, 1 Locking ring
- EXA'lence 370 HEAVY BODY
Regular and Fast Set
584g (370mL) x 1 cartridge, Mixing tip for 370mL cartridge (yellow) x 10 pieces, 1 Locking ring
- EXA Mixer 1 piece
- Dynamic Mixing Tip (Yellow) 50 pieces
- Bayonet Locking Ring 5 pieces
- Universal VPS Adhesive 1 bottle (7mL)

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3737 West 127th Street, Alsip, IL 60803, U.S.A.
GC ASIA DENTAL PTE. LTD.